





器械包 Instrument set
标准型 & 加长型 /Standard type & /Extra Long Type



苏州苏南捷迈得医疗器械有限公司
Suzhou Sunan Zimmered Medical Instrument Co.,Ltd.

地址:江苏省张家港市合兴镇昆仑山路

Add :Kunlunshan road, Hexing town, Zhangjiagang city, Jiangsu province, China

国内 China :

销售部电话 : +86-512 5850 0520
+86-186 5112 8828

国际 International :

Tel : +86-512 5691 3892
Tel/Whatsapp : +86-186 5112 8828



wechat

传真Fax: +86-512-58508580

邮编 (P.C) : 215626

邮箱E-mail:jenny.wang@zimedede.com


网址Web:www.zimedede.com

企业微信号 (Company WeChat ID) : JMDYX1990



PFNA —— 理想稳定性的最佳选择

PFNA ---Best Choice For Femur Intertrochanteric Fracture



- PFNA螺旋刀片
一个部件完成抗旋转和稳定支撑
- PFNA解剖型主钉
理想稳定 完美保证
- PFNA blade
One component achieves complete rotation
- PFNA Nail
Ideal stability perfect guarantee

苏州苏南捷迈得医疗器械有限公司
Suzhou Sunan Zimmered Medical Instrument Co.,Ltd.

产品特点 Product features



PFNA 主钉设计有
6° 外偏角
更方便从大粗隆
顶点插入
**Designed by 6°
carrying angle for
major nail**
Easier to insert



PFNA 螺旋刀 PFNA blade
一个部件完成抗旋转，稳定
支撑增加松质骨密度，恢复理想稳定。
One component can reach complete
rotation and increase the cancellous bone
density, recover soon and stable.

生物力学实验证实 Biomechanics Experiment Approved

和其他传统螺钉系统相比，PFNA 螺旋刀片和骨之间宽大的接触面积，显著提高稳定性，防止旋转和坍塌，从而避免切除（cut-out）现象。

Compared to other screw system, there's wide contact area between blade and bone, It can enhance the stability to prevent rotation and cave (cut-cut)

多个远程锁定选项 Several distal locking options

静态或动态锁定可以通过PFNA标准的小型和超小型瞄准臂进行。

Static or dynamic can be performed with the standard, small and extra small PFNA aiming arm.



静态锁定
Static locking

一步锁死恢复快
One step lock
Quick recovery



动态锁定
dynamic locking

适当避免应力遮挡
防旋转
Anti rotation
Appropriate to avoid
stress shielding effect

最佳应力分布 Optimal stress distribution

灵活的PFNA尖端可以简化插入，减轻PFNA尖端的骨骼压力。

The flexible PFNA tip eases insertion and reduces stress on the bone at the tip of the PFNA.



加长型远端特写图
Close-up view of
Lengthened distal end

PFNA —— 理想稳定性的最佳选择

PFNA --- Best Choice For Femur Intertrochanteric Fracture

产品信息 Product Information

PFNA 标准型 /Standard type

产品编号/Product No.	产品规格/Specification		
	直径/Diameter	长度/Length(mm)	颈干角/Neck-shaft Angle
4207.0918	9	170	130°
4207.0920	9	200	130°
4207.0922	9	235	130°
4207.1018	10	170	130°
4207.1020	10	200	130°
4207.1022	10	235	130°
4207.1118	11	170	130°
4207.1120	11	200	130°
4207.1122	11	235	130°

PFNA 加长型 /Lengthen type

产品编号/Product No.	产品规格/Specification		
	直径/Diameter	长度/Length(mm)	颈干角/Neck-shaft Angle
4207.0932	9	320	130° 左/右 (Left/Right)
4207.0934	9	340	130° 左/右 (Left/Right)
4207.0936	9	360	130° 左/右 (Left/Right)
4207.0938	9	380	130° 左/右 (Left/Right)
4207.0940	9	400	130° 左/右 (Left/Right)
4207.0942	9	420	130° 左/右 (Left/Right)
4207.1032	10	320	130° 左/右 (Left/Right)
4207.1034	10	340	130° 左/右 (Left/Right)
4207.1036	10	360	130° 左/右 (Left/Right)
4207.1038	10	380	130° 左/右 (Left/Right)
4207.1040	10	400	130° 左/右 (Left/Right)
4207.1042	10	420	130° 左/右 (Left/Right)
4207.1132	11	320	130° 左/右 (Left/Right)
4207.1134	11	340	130° 左/右 (Left/Right)
4207.1136	11	360	130° 左/右 (Left/Right)
4207.1138	11	380	130° 左/右 (Left/Right)
4207.1140	11	400	130° 左/右 (Left/Right)
4207.1142	11	420	130° 左/右 (Left/Right)



产品编号 product No.	长度 (mm) Length
4207. (09-11)	0/5/10

瞄准装置图例 Target device schematic

标准型远端锁定 /Stand Type distal locking



静态锁定 Static



动态锁定 Dynamic

近端锁定 /proximal locking

